# HBV surface antigen glycan isomer (HBsAgGi), correlates with HBV-DNA and **HBV-RNA** virions

K. ANGATA<sup>1</sup>, M. SOGABE<sup>1</sup>, H. NARIMATSU<sup>1</sup>, A. MURATA<sup>2</sup>, and T. GENDA<sup>2</sup> <sup>1</sup>RCMG Inc., <sup>2</sup>Juntendo University Shizuoka Hospital, Department of **Gastroenterology and Hepatology** 

### ABSTRACT:

**Background:** In HBV patients' sera, infectious particles containing HBV-DNA (Dane particles) or HBV-RNA particles are enveloped with hepatitis B virus surface antigen (HBsAg). Infectious particles contain L-, M-, and S-HBsAg, while non-infectious sub-viral particles (SVPs) mainly contain S-HBsAg. Because SVPs present much more than infectious particles, S-HBsAg is the most abundant in all HBV particles. Thus, general HBsAg test using S-HBsAg antibodies cannot distinguish infectious and non-infectious particles. Previously, we found that O-glycosylated M-HBsAg was enriched in HBV-DNA particles. In this study, we developed and examined a new antibody to detect infectious particles, which should contain O-glycosylated M-HBsAg, called as HBsAg glycan isomer (HBsAgGi). Methods and Results: Western analysis indicated that HBsAgGi antibody recognized M-HBs modified with O-glycan but not L-HBs without O-glycan on the PreS2 of genotype C. Immunoprecipitation (IP) experiments confirmed that both HBV-DNA and HBV-RNAcontaining particles were immunoprecipitated by HBsAgGi antibody. HBsAgGi localized in ER to Golgi in MHBs-expressing cells. In treatment naïve chronic hepatitis B patients, serum HBsAgGi level was significantly correlated with the HBV-DNA level (p=0.002, n=32). However, HBsAgGi level was not associated with HBV-DNA level after NA treatment, which is consistent with the IP results. **Conclusions:** O-glycosylated M-HBsAg is unique to genotype C infectious virions. New HBsAgGi antibody recognizes infectious HBV virions containing HBV-DNA or HBV-RNA, indicating that the infectious particles with HBsAgGi are generated through glycosylation pathway. Furthermore, NA treatment altered HBV-DNA and HBV-RNA status but not much HBsAgGi, suggesting that HBsAgGi is an useful monitor of infectious HBV particles during therapy.

## INTRODUCTION

Upon HBV infection, non-infectious particles are dominantly present than infectious particles containing HBV DNA (Dane particles) in patients. In contrast to non-infectious particles mainly containing S-HBsAg, Dane particles contain all L-, M-, and S-HBsAgs. General HBsAg test using antibodies for S-HBsAg detect all HBV particles. To distinguish infectious particles and non-infectious subviral particles, we generated M-HBsAg specific antibody. Because a lectin recognizing Oglycans could enrich DNA-containing particle, O-glycosylated M-HBsAg determined by mass spectrometry was targeted. This study aimed to examine a new antibody recognizing Oglycosylated M-HBsAg [HBsAg glycan isomer (HBsAgGi)] if it can detect infectious HBV particles containing HBV-DNA and/or HBV-RNA.

#### HBsAgGi antibody recognizes M-HBsAg modified with O-glycans in genotype C HBV virion



# RESULTS

- PreS2 peptides with O-glycans generated a genotype C specific antibody.
- HBsAgGi antibody recognized M-HBs but not L-HBs, which is not modified with O-glycan on the PreS2 of genotype C.
- Mutations in O-glycosylation site or removal of Oglycans resulted in decrease of the recognition by HBsAgGi antibody.
- HBsAgGi localized in ER to Golgi in M-HBs expressing cells.
- In treatment naïve CHB patients, serum HBsAgGi level was higher in HBe-positive patients (EPH) than HBenegative patients (ENH) at baseline.
- HBsAgGi levels were significantly correlated with the HBV-DNA level (n=32, *p*=0.002).
- HBsAgGi antibody can immunoprecipitate HBV-RNAcontaining particles. Percentile of HBV-RNA particles would increase after NA therapy.

HBsAgGi is specific to genotype C and present in ER to Golgi in M-HBsAg-expressing cells

CHO-K1

+ M-HBsAg

## METHODS AND MATERIALS

To generate HBsAgGi antibodies, PreS2 glycopeptides modified with O-glycans were syntehsized. To characterize the HBsAgGi antibody, western blotting, immunostaining, and ELISA were performed. To investigate clinical utility of the HBsAgGi, sera of chronic hepatitis B (CHB) patients before and after nucleos(t)ide analog (NA) treatment were analyzed by a new HBsAgGi ELISA system and compared with other HBV markers (qHBsAg, HBV-DNA, and HBcrAg). То characterize target particles of HBsAgGi, immunoprecipitated (IP) particles by HBsAgGi antibody were analyzed by reversetranscription PCR for quantifying HBV-RNA.







#### HBV particles immunoprecipitated by HBsAgGi antibody contain HBV-RNA.

HBV RNA	
---------	--

- IP (Gi) Sup

48W



#### Distribution of HBsAgGi, HBV-DNA, qHBsAg and HBcrAg at baseline in NA naïve EPH and ENH patients



#### Correlation of HBsAgGi with HBV-DNA at baseline in NA naïve CHB patients

#### HBsAgGi vs HBV DNA







#### CONCLUSIONS

- HBsAgGi antibody recognizes HBV particles in genotype specific manner dependent on the presence of Oglycans in PreS2 domain.
- HBsAgGi levels were significantly correlated with the HBV-DNA level in CHB patients during therapy.
- HBsAgGi-positive virions contain HBV-RNA.
- Thus, HBsAgGi specifically presents in infectious fraction of HBV virions containing HBV-DNA or HBV-RNA.